


Product Datasheet

SLC35B4 antibody (orb356779)

Description	SLC35B4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human UDP-xylose and UDP-N-acetylglucosamine transporter protein (174-200AA)
Target	SLC35B4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969S0
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A2780 cel...