



Product Datasheet

SLC35F5 antibody (orb356787)


Description	SLC35F5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Solute carrier family 35 member F5 protein (122-242AA)
Target	SLC35F5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WV83
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela whole cell...



Immunohistochemical staining of human lu...



Immunofluorescence analysis of Hela cell...