


Product Datasheet

SLC38A5 antibody (orb356799)

Description	SLC38A5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Sodium-coupled neutral amino acid transporter 5 protein (1-48AA)
Target	SLC38A5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WUX1
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...