

Product Datasheet

SNRNP25 antibody (Biotin) (orb356876)

Description	SNRNP25 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human U11/U12 small nuclear ribonucleoprotein 25 kDa protein (1-132AA)
Target	SNRNP25
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BV90
Expiration Date	12 months from date of receipt.