

Product Datasheet

Socs3 antibody (FITC) (orb356897)

Description	Socs3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Recombinant Rat Suppressor of cytokine signaling 3 protein (1-225AA)
Target	Socs3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	088583
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.