

Product Datasheet

SOHLH2 antibody (FITC) (orb356905)

Description	SOHLH2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Spermatogenesis- and oogenesis-specific basic helix-loop-helix-containing protein 2 protein (1-225AA)
Target	SOHLH2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NX45
Expiration Date	12 months from date of receipt.