


Product Datasheet

SPATS2L antibody (orb356915)

Description	SPATS2L antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human SPATS2-like protein (2-488AA)
Target	SPATS2L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NUQ6
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of CEM(lane 1), HeLa cells.