

## Product Datasheet

### MFSD8 antibody (HRP) (orb35693)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | MFSD8 antibody (HRP)                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                       |
| <b>Reactivity</b>          | Human                                                                                        |
| <b>Conjugation</b>         | HRP                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                        |
| <b>Immunogen</b>           | Recombinant Human Major facilitator superfamily domain-containing protein 8 protein (1-40AA) |
| <b>Target</b>              | MFSD8                                                                                        |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                               |
| <b>Form/Appearance</b>     | Liquid                                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                |
| <b>Note</b>                | For research use only                                                                        |
| <b>Isotype</b>             | IgG                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                   |
| <b>Uniprot ID</b>          | <b>Q8NHS3</b>                                                                                |
| <b>Dilution Range</b>      | ELISA                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |