


Product Datasheet

SSU72 antibody (orb356959)

Description	SSU72 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human RNA polymerase II subunit A C-terminal domain phosphatase SSU72 protein (1-194AA)
Target	SSU72
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NP77
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of PC-3 cell...