


Product Datasheet

STAC3 antibody (orb356967)

Description	STAC3 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human SH3 and cysteine-rich domain-containing protein 3 protein (1-253AA)
Target	STAC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96MF2
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat lung tissue...