


Product Datasheet

TBX15 antibody (orb357059)

Description	TBX15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human T-box transcription factor TBX15 protein (1-496AA)
Target	TBX15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96SF7
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...