

Product Datasheet

TCTEX1D2 antibody (HRP) (orb357090)

Description	TCTEX1D2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Tctex1 domain-containing protein 2 protein (1-142AA)
Target	TCTEX1D2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WW35
Expiration Date	12 months from date of receipt.