

Product Datasheet

TFEC antibody (orb357095)

Description	TFEC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Transcription factor EC protein (1-197AA)
Target	TFEC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14948
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-231 whole cell lysate.