

Product Datasheet

TIMM8B antibody (Biotin) (orb357116)

Description	TIMM8B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim8 B protein (33-83AA)
Target	TIMM8B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5J9
Expiration Date	12 months from date of receipt.