

Product Datasheet

TM4SF1 antibody (FITC) (orb357133)

Description	TM4SF1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Transmembrane 4 L6 family member 1 protein (115-161AA)
Target	TM4SF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30408
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.