



Product Datasheet

TMBIM1 antibody (orb357139)

Description	TMBIM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Protein lifeguard 3 protein (1-109AA)
Target	TMBIM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969X1
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of U251 cells...



Western blot analysis of U251 whole cel...