



Product Datasheet

TMC05A antibody (orb357143)


Description	TMC05A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Transmembrane and coiled-coil domain-containing protein 5A protein (10-192AA)
Target	TMC05A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N6Q1
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cell...



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of human te...