



Product Datasheet

TMEM192 antibody (orb357195)


Description	TMEM192 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Transmembrane protein 192 protein (193-271AA)
Target	TMEM192
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IY95
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of HeLa whole cel...



Immunoprecipitation analysis of Jurkat b...