

Product Datasheet

Tnfrsf4 antibody (orb357259)

Description	Tnfrsf4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Tumor necrosis factor receptor superfamily member 4 protein (20-211AA)
Target	Tnfrsf4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P47741
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.