

Product Datasheet

TRAPPC2L antibody (Biotin) (orb357272)

Description	TRAPPC2L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Trafficking protein particle complex subunit 2-like protein (1-139AA)
Target	TRAPPC2L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UL33
Expiration Date	12 months from date of receipt.