

Product Datasheet

TRIM51 antibody (FITC) (orb357309)

Description	TRIM51 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Tripartite motif-containing protein 51 protein (160-452AA)
Target	TRIM51
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BSJ1
Expiration Date	12 months from date of receipt.