

Product Datasheet

TRIM74 antibody (Biotin) (orb357312)

Description	TRIM74 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tripartite motif-containing protein 74 protein (47-249AA)
Target	TRIM74
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86UV6
Expiration Date	12 months from date of receipt.