

Product Datasheet

TRUB1 antibody (orb357331)

Description	TRUB1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Probable tRNA pseudouridine synthase 1 protein (2-349AA)
Target	TRUB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WWH5
Expiration Date	12 months from date of receipt.