



Product Datasheet

TTC38 antibody (orb357395)


Description	TTC38 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Tetratricopeptide repeat protein 38 protein (2-200AA)
Target	TTC38
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5R3I4
Dilution Range	WB:1:500-2000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of U251 whole cel...



Immunohistochemistry of paraffin-embedde...