

Product Datasheet

TTLL5 antibody (orb357403)

Description	TTLL5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Tubulin polyglutamylase TTLL5 protein (1-269AA)
Target	TTLL5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6EMB2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.