

Product Datasheet

UBE2O antibody (orb357435)

Description	UBE2O antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human (E3-independent) E2 ubiquitin-conjugating enzyme protein (10-277AA)
Target	UBE2O
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9C0C9
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.

