

Product Datasheet

UHRF1BP1L antibody (orb357475)

Description	UHRF1BP1L antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human UHRF1-binding protein 1-like protein (251-522AA)
Target	UHRF1BP1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0JNW5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.