



Product Datasheet

USP30 antibody (orb357495)


Description	USP30 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 30 protein (101-400AA)
Target	USP30
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q70CQ3
Dilution Range	WB:1:500-2000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CO...



WB analysis of mouse lung tissue using U...



Immunohistochemical analysis of formalin...