


Product Datasheet

UTP15 antibody (orb357503)

Description	UTP15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human U3 small nucleolar RNA-associated protein 15 homolog protein (234-518AA)
Target	UTP15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TED0
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela(lane 1), H...