

Product Datasheet

WDR18 antibody (Biotin) (orb357560)

Description	WDR18 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human WD repeat-containing protein 18 protein (231-432AA)
Target	WDR18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BV38
Expiration Date	12 months from date of receipt.