



Product Datasheet

Wnt3 antibody (orb357595)

Description	Wnt3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Proto-oncogene Wnt-3 protein (22-355AA)
Target	Wnt3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17553
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431(lane 1), H...



Western blot analysis of mouse kidney ti...