



Product Datasheet

YIPF4 antibody (orb357611)

Description	YIPF4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Protein YIPF4 protein (1-113AA)
Target	YIPF4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BSR8
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...



Western blot analysis of Jurkat whole c...