



Product Datasheet

ZDHHC6 antibody (orb357687)


Description	ZDHHC6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC6 protein (227-409AA)
Target	ZDHHC6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H6R6
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3(lane 1)...



Immunohistochemical staining of human ki...



Immunofluorescence analysis of NIH/3T3 C...