


Product Datasheet

ZFP64 antibody (orb357695)

Description	ZFP64 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Zinc finger protein 64 homolog, isoforms 1 and 2 protein (70-342AA)
Target	ZFP64
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NPA5
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver & kidney tissue.