

Product Datasheet

ZHX1-C8orf76 antibody (Biotin) (orb357704)

Description	ZHX1-C8orf76 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 protein (1-292AA)
Target	ZHX1-C8orf76
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96EF9
Expiration Date	12 months from date of receipt.