

Product Datasheet

ZNF595 antibody (HRP) (orb357770)

Description	ZNF595 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Zinc finger protein 595 protein (380-648AA)
Target	ZNF595
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IYB9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.