



Product Datasheet

ZNF614 antibody (orb357775)


Description	ZNF614 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Zinc finger protein 614 protein (1-198AA)
Target	ZNF614
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N883
Dilution Range	WB:1:500-2000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human testis tissue.



Immunohistochemical staining of human prostate tissue.



Western blot analysis of Jurkat whole cell lysate.