

Product Datasheet

ZSCAN16 antibody (FITC) (orb357817)

Description	ZSCAN16 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 16 protein (1-348AA)
Target	ZSCAN16
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H4T2
Expiration Date	12 months from date of receipt.