

Product Datasheet

Cyfra21-1 antibody (orb357828)

Description	Cyfra21-1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Cyfra21-1(BM 19.21)
Target	KRT19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08727
Expiration Date	12 months from date of receipt.