



Product Datasheet

APP antibody (orb357842)


Description	APP antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human APP Protein
Target	APP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P05067
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart & k...



Immunohistochemical staining of human te...



Immunofluorescence analysis of Hela cell...