




Product Datasheet

UBP17 antibody (orb35786)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	USP17L24
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59711
Uniprot ID	Q0WX57
NCBI	NP_001229256.1 , NP_001229255.1 , NP_001229260.1 , NP_001229259.1 , NP_001243796.1 , NP_001229261.1 , NP_001229257.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of UBP17 antibody ...



WB analysis of MAD-MB-231 cell line usin...