



Product Datasheet

Aquaporin 4 antibody (orb357880)


Description	Aquaporin 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Aquaporin-4 protein (253-323AA)
Target	AQP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55087
Dilution Range	WB:1:1000-5000, IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart & b...



Immunofluorescence analysis of U251 cell...



Immunohistochemical staining of human br...