

Product Datasheet

CD206 antibody (orb357933)

Description	CD206 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Macrophage mannose receptor 1 protein (655-1213AA)
Target	MRC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22897
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.