



Product Datasheet

NGAL antibody (orb357946)


Description	NGAL antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Neutrophil gelatinase-associated lipocalin protein
Target	LCN2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P80188
Dilution Range	IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sp...



WB analysis of NGAL transfected 293 cell...



Immunohistochemical staining of human pr...