

Product Datasheet

Rabbit Goat IgG antibody (FITC) (orb357963)

Description	Rabbit Goat IgG antibody (FITC)
Species/Host	Rabbit
Reactivity	Goat
Conjugation	FITC
Immunogen	IgG purified from goat serum
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.