

Product Datasheet

Rabbit Guinea Pig IgG antibody (orb357981)

Description	Rabbit Guinea Pig IgG antibody
Species/Host	Rabbit
Reactivity	Guinea pig
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	IgG purified from Guinea Pig serum
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.