


Product Datasheet


3X Flag antibody (orb357990)

Description	3X Flag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	DYKDDDDK(Flag) synthetic peptide conjugate to KLH
Target	Flag Tag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB:1:500-5000, IF/ICC:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.




Immunoprecipitation analysis of HEK-293

...



Western blot analysis of randomly picked...



Immunofluorescence analysis of HeLa cell...