


## Product Datasheet

### c Myc antibody (orb357992)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | c Myc antibody                                                |
| <b>Species/Host</b>        | Mouse                                                         |
| <b>Reactivity</b>          | All                                                           |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA, IF, IP, WB                                             |
| <b>Immunogen</b>           | EQKLISEEDL (Myc) synthetic peptide conjugate to KLH           |
| <b>Target</b>              | Myc tag                                                       |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG1                                                          |
| <b>Clonality</b>           | Monoclonal                                                    |
| <b>Dilution Range</b>      | WB:1:500-1:5000, IF:1:100-1:300, IP:1:1000-1:1500             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunoprecipitation analysis of HEK-293

...