


Product Datasheet

Flavivirus group antibody (orb357995)

Description	Recombinant monoclonal antibody to Flavivirus group.
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-P, WB
Immunogen	Dengue Virus type 2 antigens.
Target	Flavivirus group antigen
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.



ELISA on
Denge Virus-
Like
Particles
usin...