

Product Datasheet

Elastin antibody (orb357998)

Description

Elastin antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Elastin. Please contact us for the exact immunogen sequence. The peptide is available as orb374973.

Target

Elastin

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Application notes

This product is the same batch as product orb13391 WITHOUT BSA

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

70 kDa

Uniprot ID

[P54320](#), [P15502](#), [Q99372](#)

NCBI

[000501](#), [21](#), [31](#), [000492](#)

Entrez

[2006](#)

Dilution Range

WB: 1:200, IHC-P: 1:100, IF/ICC: 1:100

Expiration Date

12 months from date of receipt.